

Summary:

The sample is digested via Kjeldahl digestion prior to analysis in the presence of sulfuric acid, potassium sulfate, and a mercury catalyst at a final temperature of 380°C. Free ammonia and organic nitrogen compounds, such as amino acids, proteins, peptides, and other nitrogen compounds of biological origin, are converted to ammonium sulfate under these conditions. Nitrogenous compounds of some industrial wastes, such as amines, nitro compounds, hydrazones, oximes, semicarbazones, and some tertiary amines, may not be converted. The ammonium reacts with salicylate and hypochlorite in a buffered alkaline solution in the presence of sodium nitroferrocyanide (pH 12.8–13) to form the salicylic acid analog of indophenol blue. The blue-green color produced is measured at 660 nm and is proportional to the concentration of total Kjeldahl nitrogen (TKN) present in the sample.

Interferences:


Precipitation of calcium and magnesium hydroxides is eliminated by potassium sodium tartrate in the working buffer. Filter or centrifuge turbid digestates prior to analysis. Digestates with background absorbances at the analytical wavelength may interfere with the analysis.

Performance Specifications:

Range:	0.05–10 mg/L nitrogen (N)	0.20–20 mg/L N
Throughput:	51 samples/hour	51 samples/hour
Precision:		
0.05 mg/L	<8% RSD	—
0.10 mg/L	<5% RSD	—
0.20 mg/L	—	<6% RSD
1.0 mg/L	<3% RSD	<4% RSD
10 mg/L	<2% RSD	—
20 mg/L	—	<3% RSD
Method Detection Limit (MDL):	0.01 mg/L N	0.01 mg/L N

Chemicals:

Ammonium Sulfate, $(\text{NH}_4)_2\text{SO}_4$	Sodium Hydroxide, NaOH
Brij [®] -35, 30% w/v	Sodium Hypochlorite, 5.25% or 6.0% available chlorine (household bleach), NaOCl
(OI Analytical Part #A21-0110-33)	Sodium Nitroferrocyanide Dihydrate, $\text{Na}_2\text{Fe}(\text{CN})_5\text{NO}\cdot 2\text{H}_2\text{O}$
Chloroform, CHCl_3	Sodium Phosphate Dibasic Anhydrous, Na_2HPO_4
Hydrochloric Acid, concentrated, HCl	Sodium Salicylate, $\text{NaC}_7\text{H}_5\text{O}_3$
Potassium Sodium Tartrate Tetrahydrate, $\text{KNaC}_4\text{H}_4\text{O}_6\cdot 4\text{H}_2\text{O}$	Sulfuric Acid, concentrated, H_2SO_4
Potassium Sulfate, K_2SO_4	Teflon [®] or glass boiling stones
Red Mercuric Oxide, HgO	

Basic Flow Diagram:

Note: This method complies with USEPA Method 351.2.

Selected Reference: *Methods for Chemical Analysis of Water and Wastewater*; EPA-600/4-79-020; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1984; Method 351.2.

Brij is a registered trademark of ICI Americas.
Teflon is a registered trademark of E.I. du Pont de Nemours.

