


Summary:	Ammonia reacts with alkaline phenol and hypochlorite to form indophenol blue in an amount that is proportional to the ammonia concentration. The blue color is intensified with sodium nitroferrocyanide, and the absorbance is measured at 660 nm.	
Interferences:	Eliminate precipitation of calcium and magnesium hydroxides by adding sodium citrate. Filter turbid samples prior to analysis. Samples with background absorbance at the analytical wavelength may interfere.	
Performance Specifications:	Range:	0.20–250 mg/L ammonia as nitrogen (N)
	Throughput:	30 samples/hour
	Precision:	
	0.20 mg/L	<3% RSD
	250 mg/L	<2% RSD
	Method Detection Limit (MDL):	0.01 mg/L N
Chemicals:	Ammonium Sulfate, (NH ₄) ₂ SO ₄ Brij®-35, 30% w/v (OI Analytical Part #A21-0110-33) Phenol, solid or liquified, 88%, C ₆ H ₅ OH Sodium Chloride, NaCl Sodium Citrate Dihydrate, C ₆ H ₅ Na ₃ O ₇ •2H ₂ O	Sodium Hydroxide, NaOH Sodium Hypochlorite, 5.25% available chlorine (household bleach), NaOCl Sodium Nitroferrocyanide Dihydrate, Na ₂ Fe(CN) ₅ NO•2H ₂ O Sulfuric Acid, concentrated, H ₂ SO ₄

Basic Flow Diagram:



Selected References:

Methods for Chemical Analysis of Water and Wastewater; EPA-600/4-79-020; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1984; Method 351.2.

Standard Methods for the Examination of Water and Wastewater, 20th ed.; American Public Health Association: Washington, D.C., 1998.

Brij is a registered trademark of ICI Americas.