

Method Abstract**Scope**

This method is used for the determination of ammonia nitrogen in drinking water, ground water, surface water, and domestic and industrial wastes according to USEPA Method 350.1 and Standard Methods 4500-NH₃ H. This method can also be used for the determination of ammonia nitrogen in potassium chloride (KCl) extracts of soils and plants.

Summary

Ammonia reacts with alkaline phenol and hypochlorite to form indophenol blue in an amount that is proportional to the ammonia concentration. Sodium nitroferricyanide intensifies the blue color. Measure the absorbance at 640 nm. Distillation is required for regulatory compliance.

Interferences


Cyanate, which may be encountered in certain industrial effluents, hydrolyzes to some extent, even at the pH of 9.5 at which distillation is carried out. Eliminate precipitation of calcium and magnesium hydroxides by adding ethylenediaminetetraacetic acid (EDTA). Filter turbid samples prior to analysis. Samples with background absorbance at the analytical wavelength may interfere. Color intensity is sensitive to pH; standardize samples at pH 5–7 prior to analysis. Samples containing high concentrations of chlorine must be pretreated prior to analysis with sodium thiosulfate or other dechlorinating reagents.

Performance Specifications

Range:	0.01–25 mg/L
Throughput:	72 samples/hour
Precision (at 0.1 mg/L):	<1% RSD
Precision (at 1.0 mg/L):	<2% RSD
Precision (at 10 mg/L):	<5% RSD
Method Detection Limit (MDL):	0.003 mg/L
Accuracy:	97.89%

Chemicals

Ammonium Sulfate, (NH ₄) ₂ SO ₄	Sodium Hydroxide, NaOH
Brij®-35, 21% solution (part number A21-0110-33)	Sodium Hypochlorite, 5.25% available chlorine, NaOCl
Chloroform, CHCl ₃	Sodium Nitroferricyanide Dihydrate, Na ₂ Fe(CN) ₅ NO•2H ₂ O
Deionized (DI) Water, ASTM Type I or II	Sodium Thiosulfate, Na ₂ S ₂ O ₃
Ethylenediaminetetraacetic Acid, disodium salt dihydrate (EDTA), C ₁₀ H ₁₆ N ₂ Na ₂ O ₈ •2H ₂ O	Sodium Tetraborate Decahydrate, Na ₂ B ₄ O ₇ •10H ₂ O
Phenol, solid or liquefied, 88%, C ₆ H ₅ OH	Sulfuric Acid, concentrated, H ₂ SO ₄

Basic Flow Diagram

Selected Reference

Methods for the Determination of Inorganic Substances in Environmental Samples; EPA/600/R-93/100; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1993; Method 350.1.

Standard Methods for the Examination of Water and Wastewater, 21st ed.; American Public Health Association: Washington, D.C., 2005.


Figures


Figure 1. Ammonia Nitrogen Calibration (0.01–25 ppm)


Method Abstract


Figure 2. Ammonia Nitrogen Calibration Curve (0.01–25 ppm)




Figure 3. Ammonia Nitrogen Method Detection Limit (at 0.01 ppm)


Method Abstract


Figure 4. Ammonia Nitrogen Precision (at 0.10 ppm)




Figure 5. Ammonia Nitrogen ERA QC Sample Precision (ERA 1.29 ppm at 98%)


Method Abstract


Figure 6. Ammonia Nitrogen Calibration Results (0.01–25 ppm)

Method Abstract

Table 1. Ammonia Nitrogen Method Data

Parameter	Calibrant 0.01 mg/L	Calibrant 0.1 mg/L	Calibrant 1.0 mg/L	Calibrant 10.0 mg/L	ERA QC Standard 1.29 mg/L
Rep 1	0.01379	0.09543	0.09777	10.43	1.2682
Rep 2	0.01365	0.09466	0.09384	10.44	1.2578
Rep 3	0.01441	0.09624	0.09805	10.43	1.2690
Rep 4	0.01453	0.09675	0.09915	10.34	0.2593
Rep 5	0.01532	0.09618	0.09724	10.42	—
Rep 6	0.01564	0.09664	0.09788	10.45	—
Rep 7	0.01527	0.09515	0.09746	10.43	—
Rep 8	—	0.09594	0.09883	10.35	—
Rep 9	—	0.09533	0.09787	10.50	—
Rep 10	—	0.09619	0.09908	10.47	—
Average	0.01466	0.09585	0.9772	10.43	4.559
Standard Deviation	0.0007751	0.0006457	0.01436	0.04494	0.005859
% RSD	5.29	0.67	1.47	0.43	0.46
MDL	0.0024	—	—	—	—
% Accuracy	—	—	—	—	97.89%