

Method Abstract

Scope This method is used for the determination of chloride in drinking water, surface water, and domestic and industrial waste according to Standard Methods 4500-Cl⁻ E. Additionally, this method enables chloride analysis according to ISO Method 15682.

Summary Chloride reacts with mercuric thiocyanate, liberating thiocyanate ion by the formulation of soluble mercuric chloride. In the presence of ferric ion, free thiocyanate ion forms a highly colored ferric thiocyanate complex. The colored complex is measured at 480 nm.

Interferences There are no significant chemical interferences for this method. Filter turbid samples prior to analysis.


Performance Specifications

Range:	1.0–200 mg/L
Throughput:	60 samples/hour
Precision (at 10 mg/L):	<1% RSD
Precision (at 100 mg/L):	<0.5% RSD
Method Detection Limit (MDL):	0.12 mg/L
Accuracy:	98.87%

Chemicals

Brij [®] -35, 21% solution	Methanol, CH ₃ OH
Deionized (DI) Water (ASTM Type I or II)	Nitric Acid, concentrated, HNO ₃
Ferric Nitrate Nonahydrate, Fe(NO ₃) ₃ •9H ₂ O	Sodium Chloride, NaCl
Mercuric Thiocyanate, Hg(SCN) ₂	Thioacetamide, C ₂ H ₃ NS

Basic Flow Diagram



Selected Reference

Methods for Chemical Analysis of Water and Wastes; EPA-600/4-79/020; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1984; Method 325.2.

Water Quality–Determination of Chloride by Flow Analysis (CFA and FIA) and Photometric or Potentiometric Detection. *International Standard*; ISO Method 15682:2000 (E); 1st ed.; American Public Health Association: Washington, D.C., 2005.

Figures

