

## Method Abstract

**Scope** This method is used for determining phenolic materials in drinking water, surface water, and domestic and industrial wastes according to USEPA Method 420.4. Additionally, this method enables phenol index analysis according to ISO Method 14402.

**Summary** Phenol reacts with 4-aminoantipyrine and alkaline ferricyanide to form a red complex that is measured at 505 nm.

**Interferences** Interferences from sulfur compounds are eliminated by acidifying the sample to a pH of less than 4.0 with sulfuric acid, and aerating briefly by stirring and adding cupric sulfate.

Oxidizing agents such as chlorine, detected by the liberation of iodine upon acidification in the presence of potassium iodide, are removed immediately after sampling by the addition of an excess of ferrous ammonium sulfate. If chlorine is not removed, the phenolic compounds may be partially oxidized, and the results may be lower than expected.

Background contamination from reagent and sample containers can be eliminated by using glass containers. To further minimize contamination, shorten the pump tubing by cutting just before and after the Precision Pump; replace these segments of tubing with Teflon® or PET.


### Performance Specifications

Range:	5.0–500 µg/L
Throughput:	30 samples/hour
Precision (at 5.0 µg/L):	~3% RSD
Precision (at 50 µg/L):	<1% RSD
Precision (at 500 µg/L):	<1% RSD
Method Detection Limit (MDL):	0.5 µg/L
ERA QC Sample Result:	106%

### Chemicals

4-Aminoantipyrine, C <sub>11</sub> H <sub>13</sub> N <sub>3</sub> O	Phosphoric Acid, concentrated, H <sub>3</sub> PO <sub>4</sub>
Boric Acid, H <sub>3</sub> BO <sub>3</sub>	Potassium Chloride, KCl
Deionized (DI) Water, ASTM Type I or II	Potassium Ferricyanide, K <sub>3</sub> Fe(CN) <sub>6</sub>
DOWFAX® 2A1, part number A000080	Sodium Hydroxide, NaOH
Phenol, C <sub>6</sub> H <sub>5</sub> OH	

### Basic Flow Diagram



**Note** This method complies with USEPA Method 420.4.

**Selected Reference** Phenolics, Total Recoverable (Colorimetric, Semi-Automated). *Methods for the Determination of Inorganic Substances in Environmental Samples*; EPA/600/R-93/100; U.S.

Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1993; Method 420.4.

Water Quality–Determination of Phenol Index by Flow Analysis (FIA and CFA). International Standard; ISO 14402: 1999 (E); 1st ed.; Geneva, Switzerland, 1999

## Figures




Figure 1. Phenol Calibration (5.0–500  $\mu\text{g/L}$ )




Figure 2. Phenol Calibration Curve (5.0–500 µg/L)




Figure 3. Phenol Method Detection Limit (at 5 ppb)




Figure 4. Phenol Precision (at 5 ppb)




Figure 5. Phenol Precision (at 50 ppb)




Figure 6. Phenol Precision (at 500 ppb)





Figure 7. Phenol ERA QC Sample Precision (ERA 51.2 ppb)



Name	Conc	Weight
* Cal 0.00 ppb	0.000000	-0.911494
* Cal 0.00 ppb	0.000000	0.241644
* Cal 0.00 ppb	0.000000	-0.341519
* Cal 5.00 ppb	5.000000	494.140350
* Cal 5.00 ppb	5.000000	505.847015
* Cal 5.00 ppb	5.000000	476.874542
* Cal 10.0 ppb	10.000000	1117.16076
* Cal 10.0 ppb	10.000000	1141.93506
* Cal 10.0 ppb	10.000000	1108.25766
* Cal 20.0 ppb	20.000000	2013.83276
* Cal 20.0 ppb	20.000000	1953.51892
* Cal 20.0 ppb	20.000000	1954.50463
* Cal 50.0 ppb	50.000000	5075.45507
* Cal 50.0 ppb	50.000000	5081.68994
* Cal 50.0 ppb	50.000000	5031.87304
* Cal 100 ppb	100.000000	10579.0361
* Cal 100 ppb	100.000000	10531.4941
* Cal 100 ppb	100.000000	10537.8106
* Cal 200 ppb	200.000000	20925.5371
* Cal 200 ppb	200.000000	20957.8355
* Cal 200 ppb	200.000000	21037.1858
* Cal 500 ppb	500.000000	53366.0000
* Cal 500 ppb	500.000000	53329.3867
* Cal 500 ppb	500.000000	52627.0000

Calib Coef:

y=bx+a

a: (intercept) -3.8426e-01

b: 1.0483e+02

Corr Coef: 0.999952

Carryover: 0.354%

No Drift Peaks

Figure 8. Phenol Calibration Results (5.0–500 µg/L)

Table 1: Phenol Method Data

Parameter	Calibrant (5.0 µg/L)	Calibrant (5.0 µg/L)	Calibrant (50 µg/L)	Calibrant (500 µg/L)	ERA QC Standard (51.2 µg/L)
Rep 1	5.25	4.81	48.88	533.25	54.78
Rep 2	5.34	4.89	48.71	533.36	54.57
Rep 3	5.46	5.18	49.25	535.21	54.19
Rep 4	4.53	5.23	49.25	533.94	54.42
Rep 5	4.92	5.06	49.35	535.57	—
Rep 6	5.02	5.17	49.34	533.86	—
Rep 7	4.68	5.05	49.53	535.53	—
Rep 8	—	4.77	49.58	535.83	—
Rep 9	—	4.85	49.55	534.52	—
Rep 10	—	5.13	49.53	534.08	—
Average	5.028	5.055	49.188	534.515	54.489
Standard Deviation	0.3441	0.1574	0.2866	0.2930	0.2453
% RSD	6.84	3.11	0.58	0.055	0.45
MDL	1.08	—	—	—	—
% Accuracy	—	—	—	—	106.42%