


Method Abstract

Scope	This method is used for determining Total Kjeldahl Nitrogen (TKN) in drinking water, surface water, and domestic and industrial wastes according to USEPA Method 351.2 and Standard Methods 4500-N _{org} D. This method can also be used for the determination of TKN in potassium chloride (KCl) extracts of soils and plants.	
Summary	Digest the sample prior to analysis in the presence of sulfuric acid, potassium sulfate, and a mercury catalyst at a final temperature of 380 °C. Free ammonia and organic nitrogen compounds convert to ammonium sulfate under these conditions. Copper catalysts can be used, however, the green color of the catalyst interferes with the method. The ammonium reacts with salicylate and hypochlorite in a buffered alkaline solution in the presence of sodium nitroferricyanide (pH 12.8–13) to form the salicylic acid analog of indophenol blue. Measure the blue-green color produced at 660 nm.	
Interferences	Eliminate precipitation of calcium and magnesium hydroxides by adding potassium sodium tartrate in the working buffer. Filter or centrifuge turbid digestates prior to the analysis. Digestates with background absorbances at the analytical wavelength may interfere with the analysis. Mercury interference is minimized by complexation with chlorine ion.	
Performance Specifications	<div> <div>Range:</div> <div>Throughput:</div> <div>Precision (at 1.0 mg/L):</div> <div>Precision (at 5.0 mg/L):</div> <div>Method Detection Limit (MDL):</div> <div>Accuracy*:</div> </div> <div> <div>0.05–20.0 mg/L</div> <div>55 samples/hour</div> <div><1% RSD</div> <div><1% RSD</div> <div>0.01 mg/L</div> <div>99%</div> </div> <div>*Undigested; ERA (Environmental Resources Associates) WasteWatR Simple Nutrients Quality Control Sample.</div>	
Chemicals	<div>Ammonium Sulfate, (NH₄)₂SO₄</div> <div>Chloroform, CHCl₃</div> <div>Deionized Water (ASTM Type I or II)</div> <div>Hydrochloric Acid, concentrated, HCl</div> <div>Potassium Sodium Tartrate Tetrahydrate, KNaC₄H₄O₆•4H₂O</div> <div>Potassium Sulfate, K₂SO₄</div> <div>Red Mercuric Oxide, HgO</div> <div>Sodium Chloride, NaCl</div>	<div>Sodium Hydroxide, NaOH</div> <div>Sodium Hypochlorite, 5.25% available chlorine (OI Analytical does not recommend household bleach), NaOCl</div> <div>Sodium Nitroferricyanide Dihydrate, Na₂Fe(CN)₅NO•2H₂O (FW 297.95)</div> <div>Sodium Phosphate Dibasic Anhydrous, Na₂HPO₄</div> <div>Sodium Salicylate, NaC₇H₅O₃</div> <div>Sulfuric Acid, concentrated, H₂SO₄</div>

Basic Flow Diagram



Note

This method complies with USEPA Method 351.2.

Selected References

Methods for Chemical Analysis of Water and Wastewater; EPA-600/4-79-020; U.S. Environmental Protection Agency, Office of Research and Development, Environmental Monitoring and Support Laboratory: Cincinnati, OH, 1984; Method 351.2.

Standard Methods for the Examination of Water and Wastewater, 21st ed.; American Public Health Association: Washington, D.C., 2005.

Figures




Figure 1. TKN Calibration (0.05–20 ppm)




Figure 2. TKN Calibration Curve (0.05–20 ppm)




Figure 3. TKN MDL (at 0.01 ppm)




Figure 4. TKN Precision (at 1 ppm)




Figure 5. TKN Precision (at 5 ppm)




Figure 6. TKN IDAC QC Samples

DI WinFLOW (COM1--Offline) - [Analysis - finaltkn.rst]

TKN:Calibration 1: Peak 4-65

Name	Conc	Height
* cal 0.05 ppm	0.050000	760.276978
* cal 0.05 ppm	0.050000	676.993469
* cal 0.05 ppm	0.050000	697.014038
* cal 0.20 ppm	0.200000	3966.43212
* cal 0.20 ppm	0.200000	3949.56008
* cal 0.20 ppm	0.200000	4001.83886
* cal 0.50 ppm	0.500000	9246.52539
* cal 0.50 ppm	0.500000	9283.44726
* cal 0.50 ppm	0.500000	9277.96582
* cal 1.00 ppm	1.000000	17924.7539
* cal 1.00 ppm	1.000000	17916.4960
* cal 1.00 ppm	1.000000	17942.3671
* cal 2.00 ppm	2.000000	35211.9414
* cal 2.00 ppm	2.000000	35366.9339
* cal 2.00 ppm	2.000000	35409.2148
* cal 5.00 ppm	5.000000	84925.6484
* cal 5.00 ppm	5.000000	84893.3437
* cal 5.00 ppm	5.000000	84741.9296
* cal 10.00 ppm	10.000000	180415.859
* cal 10.00 ppm	10.000000	181262.296
* cal 10.00 ppm	10.000000	180561.343
* cal 20.00 ppm	20.000000	364316.629
* cal 20.00 ppm	20.000000	364135.156
* cal 20.00 ppm	20.000000	363041.531

Calib Coef:	
x=cyy+by+a	
a: (intercept)	2.0817e-03
b:	5.7112e-05
c:	-6.0600e-12
Corr Coef:	0.999916
Carryover:	0%
No Drift Peaks	

Figure 7. TKN Calibration Results (0.05–20 ppm)

Method Abstract

Table 1. Total Kjeldahl Nitrogen Precision Calculations

	1.00 mg N/L	5.00 mg N/L	0.10 mg N/L	ERA P127-739B
Rep 1	0.985	4.76	0.118	12.9
Rep 2	0.978	4.74	0.132	12.8
Rep 3	0.974	4.73	0.116	12.8
Rep 4	0.969	4.74	0.116	12.8
Rep 5	0.974	4.73	0.116	—
Rep 6	0.972	4.74	0.116	—
Rep 7	0.968	4.73	0.114	—
Rep 8	0.978	4.73	—	—
Rep 9	0.975	4.74	—	—
Rep 10	0.971	4.70	—	—
Mean	0.974	4.73	0.118	12.8
Standard Deviation	0.0005016	0.014789	0.006157	0.035151
% RSD	0.51	0.31	5.20	0.27
% Recovery	97	95	118	98
MDL	—	—	0.019	—
ERA sample recovery on non-digested, simple nutrient Ready-to-Use Quality Control Sample.				