



On-Line Distillation Cyanide by UV Digestion and Segmented Flow Analysis (SFA)

(Cartridge Part #A002599)

1.0 Scope and Application

- 1.1 This method is used for the determination of cyanide in drinking water, surface water, and domestic and industrial wastes.


Note: This method is not approved for the National Pollution Discharge Elimination System (NPDES).

- 1.2 The Method Detection Limit (MDL) of this method is 1.5 µg/L cyanide. The applicable range of the method is 5.0–500 µg/L cyanide. The range may be extended to analyze higher concentrations by sample dilution.

2.0 Summary of Method

- 2.1 Cyanide is released from cyanide complexes by UV digestion and distillation. The liberated hydrogen cyanide is converted to cyanogen chloride by reaction with chloramine-T trihydrate at a pH of less than 8. The cyanogen chloride then reacts with the pyridine-barbituric acid reagent to form a red colored complex, which is measured at 570 nm (References 15.3 and 15.5).
- 2.2 The quality of the analysis is assured through reproducible calibration and testing of the Segmented Flow Analysis (SFA) system.

- 2.3 A general flow diagram of the SFA system is shown below (see Section 17.0 for a detailed flow diagram).



3.0 Definitions

Definitions for terms used in this method are provided in Section 16.0, "Glossary of Definitions and Purposes."

4.0 Interferences

- 4.1 Interferences are eliminated or reduced prior to color formation by distillation.
- 4.2 Treat samples containing sulfide by adding cadmium nitrate powder (Section 8.6).
- 4.3 Remove oxidizing agents that decompose cyanides by adding ascorbic acid (Section 8.4).
- 4.4 Thiocyanates produce a positive interference when they are decomposed to cyanide by ultraviolet radiation (References 15.3 and 15.5).

5.0 Safety

- 5.1 The toxicity or carcinogenicity of each compound or reagent used in this method has not been fully established. Each chemical should be treated as a potential health hazard. Exposure to these chemicals should be reduced to the lowest possible level.
- 5.2 For reference purposes, a file of Material Safety Data Sheets (MSDS) for each chemical used in this method should be available to all personnel involved in this chemical analysis. The preparation of a formal safety plan is also advisable.
- 5.3 The following chemicals used in this method may be highly toxic or hazardous and should be handled with extreme caution at all times. Consult the appropriate MSDS before handling.
- 5.3.1 Barbituric Acid, $C_4H_4N_2O_3$ (FW 128.09)
 - 5.3.2 Chloramine-T Trihydrate, $CH_3C_6H_4SO_2NNaCl \cdot 3H_2O$ (FW 281.69)
 - 5.3.3 Hydrochloric Acid, concentrated, HCl (FW 36.46)
 - 5.3.4 Hypophosphorous Acid, 50% w/v, H_3PO_2 (FW 66.00)
 - 5.3.5 Phosphoric Acid, concentrated, H_3PO_4 (FW 98.00)
 - 5.3.6 Potassium Cyanide, KCN (FW 65.12)
 - 5.3.7 Potassium Phosphate Monobasic, KH_2PO_4 (FW 136.09)
 - 5.3.8 Pyridine, C_5H_5N (FW 79.10)
 - 5.3.9 Sodium Hydroxide, NaOH (FW 40.00)
 - 5.3.10 Sodium Phosphate Dibasic, Na_2HPO_4 (FW 141.96)
- 5.4 Unknown samples may be potentially hazardous and should be handled with extreme caution at all times.
- 5.5 Proper personal protective equipment (PPE) should be used when handling or working in the presence of chemicals.
- 5.6 This method does not address all safety issues associated with its use. The laboratory is responsible for maintaining a safe work environment and a current awareness file of OSHA regulations regarding the safe handling of the chemicals specified in this method.

6.0 Apparatus, Equipment, and Supplies

- 6.1 Segmented Flow Analysis (SFA) System (OI Analytical Flow Solution® IV) consisting of the following:

- 6.1.1 Model 502 Multichannel Peristaltic Pump
- 6.1.2 Random Access (RA) Autosampler
- 6.1.3 Expanded Range (ER) Photometric Detector with 5-mm path length flowcell and 660-nm optical filter
- 6.1.4 Data Acquisition System (PC or Notebook PC) with WinFLOW™ software
- 6.1.5 On-Line Distillation Cyanide Cartridge (Part #A002599)
- 6.2 Sampling equipment—Sample bottle, amber glass, with polytetrafluoroethylene (PTFE)-lined cap. Clean by washing with detergent and water, rinsing with two aliquots of reagent water, and drying by baking at 110°–150°C for a minimum of one hour.
- 6.3 Standard laboratory equipment including volumetric flasks, pipettes, syringes, etc. should all be cleaned, rinsed, and dried per bottle cleaning procedure in Section 6.2.

7.0 Reagents and Calibrants

7.1 Raw Materials

- 7.1.1 Barbituric Acid, $C_4H_4N_2O_3$ (FW 128.09)
- 7.1.2 Chloramine-T Trihydrate, $CH_3C_6H_4SO_2NNaCl \cdot 3H_2O$ (FW 281.69)
- 7.1.3 Deionized Water (ASTM Type I or II)
- 7.1.4 DOWFAX® 2A1 (Part #A000080)
- 7.1.5 Hydrochloric Acid, concentrated, HCl (FW 36.46)
- 7.1.6 Hypophosphorous Acid, 50% w/v, H_3PO_2 (FW 66.00)
- 7.1.7 Phosphoric Acid, concentrated, H_3PO_4 (FW 98.00)
- 7.1.8 Potassium Cyanide, KCN (FW 65.12)
- 7.1.9 Potassium Phosphate Monobasic, KH_2PO_4 (FW 136.09)
- 7.1.10 Pyridine, C_5H_5N (FW 79.10)
- 7.1.11 Sodium Hydroxide, NaOH (FW 40.00)
- 7.1.12 Sodium Phosphate Dibasic, Na_2HPO_4 (FW 141.96)