



## Orthophosphate in Seawater by Segmented Flow Analysis (SFA)


(Cartridge Part #A002604)

### 1.0 Scope and Application

- 1.1 This method is used for the determination of orthophosphate in seawater.
- 1.2 The Method Detection Limit (MDL) of this method is 0.009  $\mu\text{moles/L}$  orthophosphate. The applicable range of the method is 0.02–10  $\mu\text{moles/L}$  orthophosphate. The range may be extended to analyze higher concentrations by sample dilution.

### 2.0 Summary of Method

- 2.1 Orthophosphate reacts with molybdenum(VI) and antimony(III) in an acidic solution to form an antimony-phosphomolybdate complex. This complex is subsequently reduced with ascorbic acid to form a blue color, and the absorbance is measured at 880 nm (Reference 15.3).
- 2.2 The quality of the analysis is assured through reproducible calibration and testing of the Segmented Flow Analysis (SFA) system.
- 2.3 A general flow diagram of the SFA system is shown below (see Section 17.0 for a detailed flow diagram).



### 3.0 Definitions

Definitions for terms used in this method are provided in Section 16.0, "Glossary of Definitions and Purposes."

### 4.0 Interferences

- 4.1 Filter turbid samples prior to analysis.
- 4.2 The presence of more than 40 mg/L of iron(III), more than 10 mg/L of copper, or more than 10 mg/L of silica may interfere.
- 4.3 Samples with background absorbance at the analytical wavelength may interfere.

### 5.0 Safety

- 5.1 The toxicity or carcinogenicity of each compound or reagent used in this method has not been fully established. Each chemical should be treated as a potential health hazard. Exposure to these chemicals should be reduced to the lowest possible level.
- 5.2 For reference purposes, a file of Material Safety Data Sheets (MSDS) for each chemical used in this method should be available to all personnel involved in this chemical analysis. The preparation of a formal safety plan is also advisable.
- 5.3 The following chemicals used in this method may be highly toxic or hazardous and should be handled with extreme caution at all times. Consult the appropriate MSDS before handling.
  - 5.3.1 Ammonium Molybdate Tetrahydrate,  $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$  (FW 1,235.95)
  - 5.3.2 Antimony Potassium Tartrate Hemihydrate,  $\text{K}(\text{SbO})\text{C}_4\text{H}_4\text{O}_6 \cdot \frac{1}{2}\text{H}_2\text{O}$  (FW 324.92)
  - 5.3.3 Hydrochloric Acid, concentrated, HCl (FW 36.46)
  - 5.3.4 Magnesium Sulfate Heptahydrate,  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  (FW 246.48)
  - 5.3.5 Potassium Phosphate Monobasic,  $\text{KH}_2\text{PO}_4$  (FW 136.09)
  - 5.3.6 Sodium Chloride, NaCl (FW 58.44)
  - 5.3.7 Sodium Hydroxide, NaOH (FW 40.00)
  - 5.3.8 Sulfuric Acid, concentrated,  $\text{H}_2\text{SO}_4$  (FW 98.08)
- 5.4 Unknown samples may be potentially hazardous and should be handled with extreme caution at all times.

- 5.5 Proper personal protective equipment (PPE) should be used when handling or working in the presence of chemicals.
- 5.6 This method does not address all safety issues associated with its use. The laboratory is responsible for maintaining a safe work environment and a current awareness file of OSHA regulations regarding the safe handling of the chemicals specified in this method.

## **6.0 Apparatus, Equipment, and Supplies**

- 6.1 Segmented Flow Analysis (SFA) System (OI Analytical Flow Solution® IV) consisting of the following:
  - 6.1.1 Model 502 Multichannel Peristaltic Pump
  - 6.1.2 Random Access (RA) Autosampler
  - 6.1.3 Expanded Range (ER) Photometric Detector with 10-mm path length flowcell and 880-nm optical filter
  - 6.1.4 Data Acquisition System (PC or Notebook PC) with WinFLOW™ software
  - 6.1.5 Orthophosphate in Seawater Cartridge (Part #A002604)
- 6.2 Sampling equipment—Sample bottle, amber glass, with polytetrafluoroethylene (PTFE)-lined cap. Clean by washing with detergent and water, rinsing with two aliquots of reagent water, and drying by baking at 110°–150°C for a minimum of one hour.
- 6.3 Standard laboratory equipment including volumetric flasks, pipettes, syringes, etc. should all be cleaned, rinsed, and dried per bottle cleaning procedure in Section 6.2.

## **7.0 Reagents and Calibrants**

- 7.1 Raw Materials
  - 7.1.1 Ammonium Molybdate Tetrahydrate,  $(\text{NH}_4)_6\text{Mo}_7\text{O}_{24} \cdot 4\text{H}_2\text{O}$  (FW 1,235.95)
  - 7.1.2 Antimony Potassium Tartrate Hemihydrate,  $\text{K}(\text{SbO})\text{C}_4\text{H}_4\text{O}_6 \cdot \frac{1}{2}\text{H}_2\text{O}$  (FW 324.92)
  - 7.1.3 Ascorbic Acid,  $\text{C}_6\text{H}_8\text{O}_6$  (FW 176.12)
  - 7.1.4 Deionized Water (ASTM Type I or II)
  - 7.1.5 DOWFAX® 2A1 (Part #A000080)
  - 7.1.6 Hydrochloric Acid, concentrated, HCl (FW 36.46)
  - 7.1.7 Magnesium Sulfate Heptahydrate,  $\text{MgSO}_4 \cdot 7\text{H}_2\text{O}$  (FW 246.48)

7.1.8 Potassium Phosphate Monobasic,  $\text{KH}_2\text{PO}_4$  (FW 136.09)

7.1.9 Sodium Chloride,  $\text{NaCl}$  (FW 58.44)

7.1.10 Sodium Hydroxide,  $\text{NaOH}$  (FW 40.00)

7.1.11 Sulfuric Acid, concentrated,  $\text{H}_2\text{SO}_4$  (FW 98.08)