



Total Dissolved Nitrogen by Online UV/Persulfate Digestion, Segmented Flow Analysis (SFA), and Flow Injection Analysis (FIA)

(Cartridge Part #319784)


1.0 Scope and Application

- 1.1 This method is used for the determination of total dissolved nitrogen in drinking water, surface water, saline water, and domestic and industrial wastes.
- 1.2 The Method Detection Limit (MDL) of this method is 0.02 mg/L nitrogen (N). The applicable range of the method is 0.10–10 mg/L nitrogen. The range may be extended to analyze higher concentrations by sample dilution.

2.0 Summary of Method

- 2.1 Dissolved nitrogen compounds are converted to nitrate by persulfate oxidation and subsequent UV digestion in a caustic solution of boric acid (pH 14). Nitrate is reduced quantitatively to nitrite by a cadmium-copper reductor in the form of an open tube cadmium reactor. The nitrite formed is diazotized with sulfanilamide and subsequently coupled with *N*-(1-naphthyl)ethylenediamine dihydrochloride. This reaction takes place in an acidic solution (pH 2). The resulting highly colored azo dye is colorimetrically detected at 540 nm (Reference 15.4).
- 2.2 The quality of the analysis is assured through reproducible calibration and testing of the Segmented Flow Analysis (SFA) and Flow Injection Analysis (FIA) system.

- 2.3 A general flow diagram of the SFA/FIA system is shown below (see Section 17.0 for a detailed flow diagram).



3.0 Definitions

Definitions for terms used in this method are provided in Section 16.0, "Glossary of Definitions and Purposes."

4.0 Interferences

- 4.1 Eliminate interference from iron, copper, and other metals by using ethylenediaminetetraacetic acid (EDTA) in the buffer solution.
- 4.2 Norwitz and Keliher (References 15.1 and 15.2) have compiled a comprehensive study of interferences in the spectrophotometric analysis of nitrite.

5.0 Safety

- 5.1 The toxicity or carcinogenicity of each compound or reagent used in this method has not been fully established. Each chemical should be treated as a potential health hazard. Exposure to these chemicals should be reduced to the lowest possible level.
- 5.2 For reference purposes, a file of Material Safety Data Sheets (MSDS) for each chemical used in this method should be available to all personnel involved in this chemical analysis. The preparation of a formal safety plan is also advisable.
- 5.3 The following chemicals used in this method may be highly toxic or hazardous and should be handled with extreme caution at all times. Consult the appropriate MSDS before handling.
- 5.3.1 Ammonium Chloride, NH_4Cl (FW 53.50)
- 5.3.2 Ammonium Hydroxide, NH_4OH (FW 35.05)
- 5.3.3 Ammonium Sulfate, $(\text{NH}_4)_2\text{SO}_4$ (FW 132.14)

- 5.3.4 Boric Acid, H_3BO_3 (FW 61.83)
- 5.3.5 Cadmium, Cd (FW 112.40)
- 5.3.6 Chloroform, CHCl_3 (FW 119.38)
- 5.3.7 Cupric Sulfate Pentahydrate, $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ (FW 249.61)
- 5.3.8 Ethylenediaminetetraacetic Acid, Disodium Salt Dihydrate (EDTA), $\text{C}_{10}\text{H}_{16}\text{N}_2\text{Na}_2\text{O}_8 \cdot 2\text{H}_2\text{O}$ (FW 372.24)
- 5.3.9 Glycine, $\text{H}_2\text{NCH}_2\text{CO}_2\text{H}$ (FW 75.07)
- 5.3.10 Hydrochloric Acid, concentrated, HCl (FW 36.46)
- 5.3.11 *N*-(1-naphthyl)ethylenediamine Dihydrochloride, $\text{C}_{12}\text{H}_{14}\text{N}_2 \cdot 2\text{HCl}$ (FW 259.18)
- 5.3.12 Phosphoric Acid, concentrated, H_3PO_4 (FW 98.00)
- 5.3.13 Potassium Nitrate, KNO_3 (FW 101.11)
- 5.3.14 Potassium Nitrite, KNO_2 (FW 85.11)
- 5.3.15 Potassium Persulfate, $\text{K}_2\text{S}_2\text{O}_8$ (FW 270.33)
- 5.3.16 Sodium Hydroxide, NaOH (FW 40.00)
- 5.3.17 Sulfanilamide, $\text{C}_6\text{H}_8\text{N}_2\text{O}_2\text{S}$ (FW 172.21)
- 5.3.18 Urea, H_2NCONH_2 (FW 60.06)
- 5.4 Unknown samples may be potentially hazardous and should be handled with extreme caution at all times.
- 5.5 Proper personal protective equipment (PPE) should be used when handling or working in the presence of chemicals.
- 5.6 This method does not address all safety issues associated with its use. The laboratory is responsible for maintaining a safe work environment and a current awareness file of OSHA regulations regarding the safe handling of the chemicals specified in this method.

6.0 Apparatus, Equipment, and Supplies

- 6.1 Segmented Flow Analysis (SFA)/Flow Injection Analysis (FIA) System (OI Analytical Flow Solution® IV) consisting of the following:
 - 6.1.1 Model 502 Multichannel Peristaltic Pump

- 6.1.2 Random Access (RA) Autosampler
- 6.1.3 Expanded Range (ER) Photometric Detector with 5-mm path length flowcell and 540-nm optical filter
- 6.1.4 Data Acquisition System (PC or Notebook PC) with WinFLOW™ software
- 6.1.5 Total Dissolved Nitrogen Cartridge (Part #319784)
- 6.1.6 For FIA, Flow Solution IV must be equipped with the FIA option.
- 6.2 Sampling equipment—Sample bottle, amber glass, with polytetrafluoroethylene (PTFE)-lined cap. Clean by washing with detergent and water, rinsing with two aliquots of reagent water, and drying by baking at 110°–150°C for a minimum of one hour.
- 6.3 Standard laboratory equipment including volumetric flasks, pipettes, syringes, etc. should all be cleaned, rinsed, and dried per bottle cleaning procedure in Section 6.2.

7.0 Reagents and Calibrants

7.1 Raw Materials

- 7.1.1 Ammonium Chloride, NH_4Cl (FW 53.50)
- 7.1.2 Ammonium Hydroxide, NH_4OH (FW 35.05)
- 7.1.3 Ammonium Sulfate, $(\text{NH}_4)_2\text{SO}_4$ (FW 132.14)
- 7.1.4 Boric Acid, H_3BO_3 (FW 61.83)
- 7.1.5 Brij®-35, 30% w/v (Part #A21-0110-33)
- 7.1.6 Chloroform, CHCl_3 (FW 119.38)
- 7.1.7 Cupric Sulfate Pentahydrate, $\text{CuSO}_4 \cdot 5\text{H}_2\text{O}$ (FW 249.61)
- 7.1.8 Deionized Water (ASTM Type I or II)
- 7.1.9 Ethylenediaminetetraacetic Acid, Disodium Salt Dihydrate (EDTA), $\text{C}_{10}\text{H}_{16}\text{N}_2\text{Na}_2\text{O}_8 \cdot 2\text{H}_2\text{O}$ (FW 372.24)
- 7.1.10 Glycine, $\text{H}_2\text{NCH}_2\text{CO}_2\text{H}$ (FW 75.07)
- 7.1.11 Hydrochloric Acid, concentrated, HCl (FW 36.46)
- 7.1.12 N-(1-naphthyl)ethylenediamine Dihydrochloride, $\text{C}_{12}\text{H}_{14}\text{N}_2 \cdot 2\text{HCl}$ (FW 259.18)
- 7.1.13 Phosphoric Acid, concentrated, H_3PO_4 (FW 98.00)

7.1.14 Potassium Nitrate, KNO_3 (FW 101.11)

7.1.15 Potassium Nitrite, KNO_2 (FW 85.11)

7.1.16 Potassium Persulfate, $\text{K}_2\text{S}_2\text{O}_8$ (FW 270.33)

7.1.17 Sodium Hydroxide, NaOH (FW 40.00)

7.1.18 Sulfanilamide, $\text{C}_6\text{H}_8\text{N}_2\text{O}_2\text{S}$ (FW 172.21)

7.1.19 Urea, H_2NCONH_2 (FW 60.06)