




Reducing Interferences in Cyanide Analysis

William Lipps
Market Specialist – Water Analyzer Products
OI Analytical


This talk presents problems and solutions in cyanide analysis



1. What we measure




2. Problems




3. Solutions


Cyanide methods measure the various cyanide “species”




If all we had was CN^- in dilute NaOH it would be easy



Direct colorimetry does not correlate with distillation results



Direct ISE does not correlate with distilled real sample results



Cyanide methods require separation of CN from matrix


- Separated from interferences, cyanide measurement is no different than running standards.



Distillation most common technique to remove interference



Macro Distillation



**MIDI
Distillations**

Distillation actually creates CN interferences

- **Boiling acid**
- **Automated UV-Distillation**
 - Boiling acid





Interferences with distillation are in almost every sample

- **Thiocyanate**
- **Sulfur**
- **Thiosulfate/Sulfite**
- **Oxidizers**

The left side of the slide features a vertical blue bar with a white chromatogram line. Along this line are four boxes containing chemical symbols: CN^- , NH_3 , PO_4 , and NO_3^- .

Thiocyanate + Nitrate results in positive bias

- The addition of Sulfamic acid does not sufficiently reduce this interference.
 - A real POTW sample with 0.1 mg/L SCN^- and 63.5 mg/L NO_3^- detected total CN^- at **0.10 mg/L** even after the addition of Sulfamic Acid

Sulfur compounds react rapidly with CN

- **Elemental Sulfur**
 - $8\text{CN}^- + \text{S}_8 \rightarrow \text{SCN}^-$
- **Metal Sulfides**
 - Cu_2S , FeS , PbS , CuFeS_2 , CdS , ZnS , etc.
 - S reacts with CN^- to form SCN^-



Thiosulfate reacts with cyanide during distillation

- **0.200 mg/L CN^- + 200 mg/L $\text{S}_2\text{O}_3^{2-}$**
 - Cyanide Found = 0.160 mg/L
 - Recovery = 80%*
- * Double Chloramine T added, or recovery would be lower.

Sulfite reacts rapidly with CN in basic solutions

- **0.200 mg/L CN⁻ + 200 mg/L SO₃⁻²**
 - Cyanide Found = 0.000 mg/L
 - Recovery = 0%
- This reaction occurs in absorber solution, or in preserved sample



There is no way to “know” if sulfur compounds are present

- **No “spot” tests that adequately detect the sulfur compounds**
- **Sodium sulfite and sodium thiosulfate are both added to samples for dechlorination.**

CN

NH₃

PO₄

NO₃



Oxidizers destroy cyanide before or during distillation

- Hypochlorite
- Peroxide
- Caro's Acid
- Chloramines

The left side of the slide features a vertical blue bar with a white chromatogram line. Along this line are four boxes containing chemical formulas: CN, NH₃, PO₄, and NO₃.

Footnote 6 (MUR 2007) allows other methods to be used

- **More accurate?**
 - Spikes?
- **“challenge matrix” distilled**

Matrix spikes cannot be used to demonstrate accuracy

Method	Amount Detected (ppb)	Recovery
335.4	32	98 %
335.3	16	98 %

Both methods detected CN^- in a synthetic sample with no CN^- .



Verify Accuracy Using Interference Free Methods

- Use methods demonstrated by literature and multiple users to be interference free
 - OIA 1677 or ASTM D6888-04
 - ASTM D 7284-08
 - OIA 1678 (ASTM D7511-09)

CN

NH₃


PO₄

NO₃

The old versus the new




**Purging and
boiling in acid**



**Diffusion at
room
temperature**

Electrochemistry techniques integrate matrix removal

- Very sensitive with large dynamic range.





Unlike colorimetry, GD amperometry is easy to visualize

- $\text{CN}^- + \text{H}^+ \rightarrow \text{HCN}$
- $\text{HCN} + \text{OH}^- \rightarrow \text{CN}^- + \text{H}_2\text{O}$
- $\text{Ag} + 2\text{CN}^- \rightarrow \text{Ag}(\text{CN})_2^- + \text{e}^-$



measure


CN

NH₃

PO₄

NO₃

This flow diagram illustrates the simplicity of GD-amperometry





A series of GD-amperometry CN methods to meet needs


- **Free cyanide**
- **Available cyanide**
- **Total distilled cyanide**
- **Total non-distilled cyanide**




The only fully automated free cyanide method

Method	Description	Measurement
ASTM D 7237	FIA	Gas Diffusion-Amperometry

Direct colorimetry is not measuring free cyanide




The free cyanide method by GD-amperometry



ASTM D7237-06

Real Data Comparison – Free CN



GD-amperometry methods for available cyanide – no CATC!


Descriptive Name	Method Number	Description	Measurement
Available Cyanide	OIA 1677	Ligand Exchange / Flow Injection Analysis	Gas Diffusion - Amperometry
	ASTM D 6888	Ligand Exchange / Flow Injection Analysis	Gas Diffusion - Amperometry



No distillation or pyridine required


Sulfide does not interfere

Sulfide Interference on Cyanide Methods
Recovery of 0.1 ppm Cyanide in Samples Containing Sulfide




Quantitative recovery compared to other methods such as CATC

Comparison of Recoveries by Available Cyanide Methods
Concentration Dependent Recovery by Cyanide Species



Precise and more accurate than WAD methods

Comparison OIA 1677 vs WAD and WAD-Kelada Distillations
Real World Samples Fortified at 100 ppb CN



GD-amperometry available CN has fewer interferences

CATC	WAD	OIA 1677
N-organics	Excessive Iron Cyanide	None
SCN, NH ₃ , NO ₂	Concentration Dependent	
S ₂ O ₃ , H ₂ O ₂		
Concentration Dependent		

GD-amperometry available CN saves time and labor

	CATC	WAD	OIA 1677
Sample Preparation	2 distillations 2 – 3 hours	1 distillation 2 – 3 hours	No distillation
Analysis	1 – 2 minutes	1 – 2 minutes	1 – 2 minutes
Total Time	3 – 4 hours	3 – 4 hours	1 – 2 minutes



Less manipulation means better Available CN data

- **No distillation**
- **0.5 ppb MDL**
- **Up to 90 samples per hour**
- **Ease of Operation**
- **Very simple chemistry**

CN

NH₃


PO₄

NO₃

Easy to use, understand, and environmentally friendly



No pyridine



**OIA 1677, ASTM D6888-04 or
ASTM D7284-08**



Total cyanide methods using manual distillation

Descriptive Name	Method Number	Description	Measurement
Total Cyanide	EPA 335.4	Midi Distillation – MgCl_2	Automated Colorimetry
	ASTM D 7284	Midi / Micro Distillation – MgCl_2	Gas Diffusion - Amperometry

CN

NH_3

PO_4

NO_3



Most total cyanide analyses are by EPA 335.4 or similar

- Manual distillation
- Prolonged heating (125 °C) , strong acid (pH <2)
- Purging into basic absorber solution
- Colorimetry

CN

NH₃

PO₄

NO₃



ASTM D7284 manual distillation GD-amperometry method

- **Distill samples**
- **Use GD-amperometry**
- **No pyridine**
- **Fewer interferences**

CN


NH₃

PO₄

NO₃

Once again, sulfide does not interfere with GD-amperometry

Sulfide Interference on Cyanide Methods
Recovery of 0.1 ppm Cyanide in Samples Containing Sulfide



**But wait, I thought
distillation was bad?**




Comparison of Measurements with Interferences Present

Sample (Total CN)	EPA 335.4	ASTM D7284
20 ppb	18.7	18.4
20 ppb +SCN + NO ₃	227	229
200 ppb + SO ₃	147	159
200 ppb + Ascorbate	152	154
200 ppb + Ascorbate + Ammonia	73	71
200 ppb + Ascorbate + Ammonia + OCL ⁻	47	49

What exactly do “total” cyanide methods measure?

- **Iron Cyanide + Available CN**



Iron cyanide is not toxic.

- Sunlight causes iron cyanide to release HCN
- Sunlight = UV irradiation



Automated total cyanide methods use UV to liberate HCN from Fe


Descriptive Name	Method Number	Description	Measurement
Total Cyanide	ASTM D4374 (Kelada 01)	High power UV- Auto distillation Alkaline pH	Automated colorimetry
	EPA 335.3	Low power UV- Auto distillation pH <2	Automated Colorimetry
	OIA 1678/ASTM D7511	Low power UV- pH <2	Gas Diffusion - Amperometry

Comparison of Kelada and ASTM D7511-09

	Kelada 01	ASTM D7511
Pump Tubes	15	5
Reagents	Pyridine	No Pyridine
Distillation	Yes	No
SCN Interaction	0.25 – 0.5 %	0.01 – 0.03 %


For total cyanide, simply add UV irradiation to the GD method

No pyridine




OIA1678 / ASTM D7511-09

Same results with and without flash distillation




Same, but more precise results as manual distillation

Comparison of OIA 1678 and EPA 335.4
Real World Samples



ASTM D7511 and D7284 get the same result if no interferences

Comparison Distillation with Non Distillation Methods
Real World Samples



Comparison of Total CN methods

	335.4	ASTM D7284	ASTM D7511
Sample Preparation	2 – 3 hour distillation	1 – 3 hour distillation	No distillation
Analysis	1 – 2 minutes	1 – 2 minutes	1 – 2 minutes
Total Time	3 – 4 hours	2 – 4 hours	1 – 2 minutes

In summary, distillation/colorimetry should be replaced

